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Ziegelei 1

D-72336 Balingen

E-Mail: [info@kern-sohn.com](mailto:info@kern-sohn.com)

Tel: +49-[0]7433- 9933-0

Fax: +49-[0]7433-9933-149

Internet: [www.kern-sohn.com](http://www.kern-sohn.com)

# **Zusatzbeschreibung Supplementary instructions Instructions supplémentaires**

## **KERN KFE-TM**

Version 1.0  
10/2011

Deutsch

English

Français

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---

## 1 Adjustment

As the acceleration value due to gravity is not the same at every location on earth, each display unit with connected weighing plate must be coordinated - in compliance with the underlying physical weighing principle - to the existing acceleration due to gravity at its place of location (only if the weighing system has not already been adjusted to the location in the factory). This adjustment process must be carried out for the first commissioning, after each change of location as well as in case of fluctuating environment temperature. To receive accurate measuring values it is also recommended to adjust the display unit periodically in weighing operation.

- The adjustment is locked for verified balances.

In order to unlock the access, the seal must be destroyed and the jumper on the printed circuit board must be fitted.

Attention:

After destruction of the seal the weighing system must be re-verified by an authorised agency and a new verification wire/seal mark fitted before it can be reused for applications subject to verification.

- The weight to be used depends on the capacity of the scale. Carry out adjustment as near as possible to the scale's maximum weight. Info about test weights can be found on the Internet at: <http://www.kern-sohn.com>.
- Observe stable environmental conditions. Stabilisation requires a certain warm-up time.

## Call-up menu

|                                                                                                                                                                                                                                                                                                                                                               |                                                         |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------|
| ⇒ In weighing mode press  and  at the same time and the first menu block <i>F0 H-L</i> will appear.                                                                                         | <i>F0 H-L</i>                                           |
| ⇒ Press  repeatedly until <i>ProG</i> is displayed.                                                                                                                                                                                                                          | <i>ProG</i>                                             |
| ⇒ Press  , password query <i>P1n</i> will appear.                                                                                                                                                                                                                            | <i>P1n</i>                                              |
| ⇒ Press  ,  ,  subsequently until the first menu item <i>P1rEF</i> is shown.                               | <i>P1rEF</i>                                            |
| ⇒ Press  repeatedly until <i>P2n0d</i> is displayed.                                                                                                                                                                                                                         | <i>P2n0d</i>                                            |
| ⇒ Confirm by  .<br>⇒ Press  repeatedly until the currently balance typ will be displayed.<br><br><i>S1Gr</i> = Singel range<br><i>dUA rA</i> = Dual range<br><i>dUA in</i> = Multi interval | <i>S1Gr</i><br>⇕<br><i>dUA rA</i><br>⇕<br><i>dUA in</i> |
| ⇒ Confirm by  .                                                                                                                                                                                                                                                            | <i>dec.</i>                                             |
| ⇒ Press  repeatedly until <i>CAL</i> will be displayed.                                                                                                                                                                                                                    | <i>CAL</i>                                              |
| ⇒ Confirm by  and select by  .<br><i>nonLin</i> = Justierung<br><i>LinEAr</i> = Linearisierung                                                                                          | <i>nonLin</i><br>⇕<br><i>LinEAr</i>                     |

## Procedure

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |                                                                                                                                                                                 |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>⇒ Confirm menu setting <i>nonLin</i> by . Ensure that there are no objects on the weighing plate.</p> <p>⇒ Wait for stability display, then press .</p>                                                                                                                                                                 | <br>↓<br> |
| <p>⇒ The currently set adjustment weight will be displayed.</p>                                                                                                                                                                                                                                                                                                                                                                                                                              |                                                                                              |
| <p>⇒ Either use the displayed adjustment weight or change it with ,  and , the active digit is flashing.</p> <p>⇒ Confirm by , „LoAd“ will be shown.</p> |                                                                                              |
| <p>⇒ Carefully place adjusting weight in the centre of the weighing plate. Wait for stability display, then press .</p>                                                                                                                                                                                                                                                                                     |                                                                                              |
| <p>⇒ After the adjustment the balance will carry out a self-test. Remove adjusting weight <b>during</b> selftest, balance will return into weighing mode automatically. An adjusting error or incorrect adjusting weight will be indicated by the error message; repeat adjustment procedure.</p>                                                                                                                                                                                            |                                                                                            |

## 2 Linearization

Linearity shows the greatest deviation of a weight display on the scale to the value of the respective test weight according to plus and minus over the entire weighing range.

If linearity deviation is discovered during a testing instrument control, you can improve this by means of linearization.

- Carrying out linearization is restricted to specialist staff possessing well acquainted with the workings of weighing scales.
- The linearisation is locked for verified weighing systems. In order to unlock the access, the seal must be destroyed and the jumper on the printed circuit board must be fitted.  
Attention:  
After destruction of the seal the weighing system must be re-verified by an authorised agency and a new verification wire/seal mark fitted before it can be reused for applications subject to verification.
- The test weights to be used must be adapted to the weighing scale's specifications.
- Observe stable environmental conditions. Stabilisation requires a certain warm-up time.
- After successful linearization you will have to carry out calibration; see chapter "testing instruments control"

|                                                                                                                                                                                                                                                                                                                |                                                                                                                                                                                                                                                                   |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p><b>Procedure:</b></p> <p>⇒ Call-up menu setting <i>LinEAr</i>, see chap. 1.</p> <p>⇒ Acknowledge with . Ensure that there are no objects on the weighing plate.</p>                                                        | <br><br> |
| <p>⇒ Wait for stability display „STABLE“, then press . When “LoAd 1” is displayed, put the first adjustment weight (1/3 max) carefully in the centre of the weighing platform.</p>                                            |                                                                                                                                                                                |
| <p>⇒ Wait for stability display „STABLE“, then press . When “LoAd 2” is displayed, put the second adjustment weight (2/3 max) carefully in the centre of the weighing platform.</p>                                           |                                                                                                                                                                                |
| <p>⇒ Wait for stability display „STABLE“, then press . When “LoAd 3” is displayed, put the third adjustment weight (1/3 max) carefully in the centre of the weighing platform.</p>                                            |                                                                                                                                                                                |
| <p>⇒ Wait for stability display „STABLE“, then press .</p>                                                                                                                                                                  |                                                                                                                                                                              |
| <p>⇒ After the adjustment the balance will carry out a self-test. Remove adjusting weight <b>during</b> selftest, the appliance will return into weighing mode automatically.</p> <p>An adjusting error or incorrect adjusting weight will be indicated by the error message; repeat adjustment procedure.</p> |                                                                                                                                                                              |

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